

Product Name: CLPP Rabbit Monoclonal antibody
Catalog #: AMRe01834



Summary

Production Name	CLPP Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CLPP
Alternative Names	Endopeptidase Clp
Gene ID	8192
SwissProt ID	Q16740.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 30 kDa; Observed MW: 30 kDa

Background

Clp cleaves peptides in various proteins in a process that requires ATP hydrolysis. Clp may be responsible for a fairly

Product Name: CLPP Rabbit Monoclonal antibody
Catalog #: AMRe01834




general and central housekeeping function rather than for the degradation of specific substrates.

Research Area

Signal Transduction

Image Data



Western blot analysis of CLPP in K562, rat Brain, 3T3, HeLa lysates using CLPP antibody.

Note

For research use only.