

Product Name: Cullin 1 Rabbit Monoclonal antibody
Catalog #: AMRe01861



Summary

Production Name	Cullin 1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	CUL1
Alternative Names	CUL1; Cullin1
Gene ID	8454
SwissProt ID	Q13616.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 90 kDa; Observed MW: 90 kDa

Background

Product Name: Cullin 1 Rabbit Monoclonal antibody
Catalog #: AMRe01861




Core component of multiple cullin-RING-based SCF (SKP1-CUL1-F-box protein) E3 ubiquitin-protein ligase complexes, which mediate the ubiquitination of proteins involved in cell cycle progression, signal transduction and transcription.


Research Area

Cell Biology


Image Data



Western blot analysis of Cullin 1/CUL1 in MEF, 3T3 lysates using Cullin 1 antibody.



Western blot analysis of Cullin 1/CUL1 in K562, C6 lysates using Cullin 1/CUL1 antibody.



Immunocytochemistry analysis of Cullin 1(green) in HeLa using Cullin 1 antibody, and DAPI(blue)

Product Name: Cullin 1 Rabbit Monoclonal antibody
Catalog #: AMRe01861



Note

For research use only.