

**Product Name: Cytochrome C Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01881**

---

## Summary

<b>Production Name</b>	Cytochrome C Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB, ICC/IF, IP
<b>Reactivity</b>	Human, Mouse, Rat

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	CYCS
<b>Alternative Names</b>	CYCS; CYC; Cytochrome c
<b>Gene ID</b>	54205
<b>SwissProt ID</b>	P99999.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
<b>Molecular Weight</b>	Calculated MW: 12 kDa; Observed MW: 12 kDa

**Product Name: Cytochrome C Rabbit Monoclonal Antibody**  
**Catalog #: AMRe01881**

---


## Background

CYCS Electron carrier protein. The oxidized form of the cytochrome c heme group can accept an electron from the heme group of the cytochrome c1 subunit of cytochrome reductase. Cytochrome c then transfers this electron to the cytochrome oxidase complex, the final protein carrier in the mitochondrial electron-transport chain.


## Research Area

Cardiovascular

## Image Data



Immunocytochemistry analysis of Cytochrome C (green) in K562 using Cytochrome C antibody, and DAPI (blue).



Western blot analysis of Cytochrome C in K562, rat Brain, C6, 3T3, HeLa lysates using Cytochrome C antibody.

## Note

For research use only.