

Product Name: DAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe01893



Summary

Production Name	DAP1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	DAP
Alternative Names	Death-associated protein 1; DAP-1
Gene ID	1611
SwissProt ID	P51397.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 11 kDa; Observed MW: 11 kDa

Background

Negative regulator of autophagy. Involved in mediating interferon-gamma-induced cell death.


Product Name: DAP1 Rabbit Monoclonal Antibody
Catalog #: AMRe01893



Research Area

Cell Biology

Image Data



Western blot analysis of DAP1 in MCF-7, C6, 3T3 lysates using DAP1 antibody.

Note

For research use only.