

Product Name: Destrin Rabbit Monoclonal Antibody
Catalog #: AMRe01904



Summary

Production Name	Destrin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	DSTN
Alternative Names	ACTDP; ADF; corn1; Destrin; DSN; Dstn; HEL32; Sid 23; sid23p
Gene ID	11034
SwissProt ID	P60981.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 19 kDa; Observed MW: 19 kDa

Background

Product Name: Destrin Rabbit Monoclonal Antibody
Catalog #: AMRe01904



Actin-depolymerizing protein. Severs actin filaments (F-actin) and binds to actin monomers (G-actin). Acts in a pH-independent manner.


Research Area

Cell Biology

Image Data



Western blot analysis of Destrin in C6 lysates using Destrin antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Destrin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.