

Product Name: DNAJA4 Rabbit Monoclonal antibody
Catalog #: AMRe01917



Summary

Production Name	DNAJA4 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	DNAJA4
Alternative Names	MST104; MSTP104; PRO1472
Gene ID	55466
SwissProt ID	Q8WW22.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 45 kDa; Observed MW: 45 kDa

Background

Product Name: DNAJA4 Rabbit Monoclonal antibody
Catalog #: AMRe01917




Enables chaperone binding activity and unfolded protein binding activity. Involved in several processes, including negative regulation of endothelial cell migration; negative regulation of inclusion body assembly; and protein refolding. Located in cytosol and membrane.


Research Area

Signal Transduction

Image Data



Western blot analysis of DNAJA4 in HeLa lysates using DNAJA4 antibody.



Western blot analysis of DNAJA4 in C6, rat Brain lysates using DNAJA4 antibody.

Note

For research use only.