

Product Name: DPS1 Rabbit Monoclonal Antibody
Catalog #: AMRe01925



Summary

Production Name	DPS1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	PDSS1
Alternative Names	Trans-prenyltransferase 1; TPT 1; Solanesyl-diphosphate synthase subunit 1; Decaprenyl-diphosphate synthase subunit 1
Gene ID	23590
SwissProt ID	Q5T2R2.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 46 kDa; Observed MW: 46 kDa

Background

Product Name: DPS1 Rabbit Monoclonal Antibody
Catalog #: AMRe01925




Decaprenyl-diphosphate synthase subunit 1

Research Area

Signal Transduction

Image Data



Western blot analysis of DPS1 in K562 lysates using DPS1 antibody.

Note

For research use only.