

Product Name: FKBP51 Rabbit Monoclonal Antibody
Catalog #: AMRe01988



Summary

Production Name	FKBP51 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	FKBP5
Alternative Names	P54; AIG6; FKBP51; FKBP54; PPlase; Ptg-10
Gene ID	2289
SwissProt ID	Q13451.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 51 kDa; Observed MW: 51 kDa

Background

Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.


Product Name: FKBP51 Rabbit Monoclonal Antibody
Catalog #: AMRe01988



Research Area

Immunology

Image Data



Western blot analysis of FKBP51 in CHO-K1 lysates using FKBP51 antibody.

Note

For research use only.