

Product Name: FRA2 Rabbit Monoclonal Antibody
Catalog #: AMRe02002



Summary

Production Name	FRA2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	FOSL2
Alternative Names	FOSL2; FRA2; Fos-related antigen 2; FRA-2
Gene ID	2355
SwissProt ID	P15408.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35-45 kDa

Background

Controls osteoclast survival and size (By similarity). As a dimer with JUN, activates LIF transcription (By similarity). Activates

Product Name: FRA2 Rabbit Monoclonal Antibody
Catalog #: AMRe02002




CEBPB transcription in PGE2-activated osteoblasts.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of FRA2 in K562, HeLa lysates using FRA2 antibody.

Note

For research use only.