

Summary

Production Name	FUBP1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	FUBP1
Alternative Names	DNA helicase V; FBP; FUBP; Fubp1; Fubp4; FUSE-binding protein 1; HDH V
Gene ID	8880
SwissProt ID	Q96AE4.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 68 kDa; Observed MW: 74 kDa

Background

Product Name: FUBP1 Rabbit Monoclonal antibody
Catalog #: AMRe02008




Regulates MYC expression by binding to a single-stranded far-upstream element upstream of the MYC promoter. May act both as activator and repressor of transcription.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Immunocytochemistry analysis of FUBP1 (green) in 293 using FUBP1 antibody, and DAPI (blue).




Western blot analysis of FUBP1/FBP in mouse spleen lysates using FUBP1/FBP antibody.



Western blot analysis of FUBP1/FBP in K562, rat Brain, HeLa lysates using FUBP1/FBP antibody

Product Name: FUBP1 Rabbit Monoclonal antibody
Catalog #: AMRe02008



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using FUBP1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.