

Summary

Production Name	GSTK1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	GSTK1
Alternative Names	GST; GST class kappa; GST13; Gstk1; hGSTK1
Gene ID	373156
SwissProt ID	Q9Y2Q3.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 25 kDa; Observed MW: 25 kDa

Background

Product Name: GSTK1 Rabbit Monoclonal Antibody
Catalog #: AMRe02044




Significant glutathione conjugating activity is found only with the model substrate, 1-chloro-2,4-dinitrobenzene (CDNB).


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of GSTK1 in HeLa lysates using GSTK1 antibody.



Immunocytochemistry analysis of Glutathione S Transferase kappa 1(green) in HeLa using Glutathione S Transferase kappa 1 antibody, and DAPI(blue)

Note

For research use only.