

Product Name: GST3 Rabbit Monoclonal antibody
Catalog #: AMRe02065



Summary

Production Name	GST3 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	GSTP1
Alternative Names	GSTP1; FAES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1
Gene ID	2950
SwissProt ID	P09211.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Product Name: GST3 Rabbit Monoclonal antibody
Catalog #: AMRe02065




Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.


Research Area

Cardiovascular


Image Data



Immunocytochemistry analysis of GST3 (green) in A549 using GST3 antibody, and DAPI (blue).



Western blot analysis of GST3/GST pi in K562, C6, 3T3, HeLa lysates using GST3/GST pi antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer tissue using GST3/GST pi antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: GST3 Rabbit Monoclonal antibody
Catalog #: AMRe02065



Note

For research use only.