

Product Name: HARS Rabbit Monoclonal antibody
Catalog #: AMRe02070



Summary

Production Name	HARS Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	HARS1
Alternative Names	HRS; CMT2W; USH3B
Gene ID	3035
SwissProt ID	P12081.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 57 kDa; Observed MW: 57 kDa

Background

Cytoplasmic histidine--tRNA ligase (Probable). Plays a role in axon guidance.


Product Name: HARS Rabbit Monoclonal antibody
Catalog #: AMRe02070



Research Area

Immunology

Image Data



Western blot analysis of HARS in K562, HeLa lysates using HARS antibody.

Note

For research use only.