

Product Name: Histone H1.0 Rabbit Monoclonal antibody
Catalog #: AMRe02082

Summary

Production Name	Histone H1.0 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	H1-0
Alternative Names	Histone H1.0; Histone H1(0); Histone H1.0; N-terminally processed; H1F0; H1FV; Histone H5
Gene ID	3005
SwissProt ID	P07305.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 21 kDa; Observed MW: 30 kDa

Product Name: Histone H1.0 Rabbit Monoclonal antibody
Catalog #: AMRe02082


Background

Histone H1.0 is a lysine rich member of the H1 family of linker histones. The H1 family of proteins interacts with linker DNA between nucleosomes and mediates compaction into higher order chromatin. Histones H1 are necessary for the condensation of nucleosome chains into higher-order structures. The H1F0 histones are found in cells that are in terminal stages of differentiation or that have low rates of cell division.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Histone H1.0 in C6, CHO-K1 lysates using Histone H1.0 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Histone H1.0 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.