

Product Name: HLTF Rabbit Monoclonal Antibody
Catalog #: AMRe02096



Summary

Production Name	HLTF Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	HLTF
Alternative Names	HIP116; HIP116A; HLTF 1; Hltf; HLTF1; p113; RNF80; SMARCA3; SMARCA3; SMARCA3;
Gene ID	6596
SwissProt ID	Q14527.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 114 kDa; Observed MW: 114 kDa

Product Name: HLTF Rabbit Monoclonal Antibody
Catalog #: AMRe02096




Background

Has both helicase and E3 ubiquitin ligase activities. Possesses intrinsic ATP-dependent nucleosome-remodeling activity; This activity may be required for transcriptional activation or repression of specific target promoters (By similarity).


Research Area

Cell Biology

Image Data



Western blot analysis of HLTF in HeLa lysates using HLTF antibody.



Immunocytochemistry analysis of HLTF(green) in HeLa using HLTF antibody, and DAPI(blue)

Note

For research use only.