

Product Name: IQGAP2 Rabbit Monoclonal Antibody
Catalog #: AMRe02168



Summary

Production Name	IQGAP2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	IQGAP2
Alternative Names	Ras GTPase-activating-like protein IQGAP2
Gene ID	10788
SwissProt ID	Q13576.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 181 kDa; Observed MW: 181 kDa

Background

Binds to activated CDC42 and RAC1 but does not seem to stimulate their GTPase activity. Associates with calmodulin.


Product Name: IQGAP2 Rabbit Monoclonal Antibody
Catalog #: AMRe02168



Research Area

Signal Transduction

Image Data



Western blot analysis of IQGAP2 in Ramos lysates using IQGAP2 antibody.

Note

For research use only.