

**Product Name: Lactate Dehydrogenase C Rabbit  
Monoclonal Antibody  
Catalog #: AMRe02199**

---

## Summary

<b>Production Name</b>	Lactate Dehydrogenase C Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB
<b>Reactivity</b>	Human

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	LDHC
<b>Alternative Names</b>	CT32; LDH3; LDHX
<b>Gene ID</b>	3948
<b>SwissProt ID</b>	P07864.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000
<b>Molecular Weight</b>	Calculated MW: 36 kDa; Observed MW: 36 kDa

## Background

---

**Product Name: Lactate Dehydrogenase C Rabbit  
Monoclonal Antibody  
Catalog #: AMRe02199**

---




Possible role in sperm motility.

## Research Area

Signal Transduction

## Image Data



Western blot analysis of Lactate Dehydrogenase C in HeLa lysates using Lactate Dehydrogenase C antibody.

## Note

For research use only.