

Product Name: N4BP1 Rabbit Monoclonal antibody
Catalog #: AMRe02304



Summary

Production Name	N4BP1 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	N4Bp1
Alternative Names	C81621; A1481586
Gene ID	80750.0
SwissProt ID	Q6A037.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 99 kDa; Observed MW: 99 kDa

Background

Potent suppressor of cytokine production that acts as a regulator of innate immune signaling and inflammation.


Product Name: N4BP1 Rabbit Monoclonal antibody
Catalog #: AMRe02304



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of N4BP1 in K562, HeLa lysates using N4BP1 antibody.

Note

For research use only.