

Product Name: PFKL Rabbit Monoclonal Antibody
Catalog #: AMRe02428



Summary

Production Name	PFKL Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	PFKL
Alternative Names	6-phosphofructokinase type B; PFK-B; Phosphohexokinase
Gene ID	5211
SwissProt ID	P17858.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 85 kDa; Observed MW: 85 kDa

Background

Catalyzes the phosphorylation of D-fructose 6-phosphate to fructose 1,6-bisphosphate by ATP, the first committing step of

Product Name: PFKL Rabbit Monoclonal Antibody
Catalog #: AMRe02428




glycolysis.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of PFKL in C6, 3T3 lysates using PFKL antibody.

Note

For research use only.