

Summary

| | |
|------------------------|---------------------------------|
| Production Name | PON2 Rabbit Monoclonal Antibody |
| Description | Rabbit Monoclonal antibody |
| Host | Rabbit |
| Application | WB,IHC-P,IP |
| Reactivity | Human |

Performance

| | |
|---------------------|--|
| Conjugation | Unconjugated |
| Modification | Unmodified |
| Isotype | IgG |
| Clonality | Monoclonal |
| Form | Liquid |
| Storage | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| Buffer | 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA |
| Purification | Affinity Purification |

Immunogen

| | |
|--------------------------|-----------------------------------|
| Gene Name | PON2 |
| Alternative Names | A esterase 2; Paraoxonase 2; PON2 |
| Gene ID | 5445 |
| SwissProt ID | Q15165. |

Application

| | |
|-------------------------|---|
| Dilution Ratio | WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20 |
| Molecular Weight | Calculated MW: 39 kDa; Observed MW: 40-42 kDa |

Background

Product Name: PON2 Rabbit Monoclonal Antibody
Catalog #: AMRe02459




Capable of hydrolyzing lactones and a number of aromatic carboxylic acid esters. Has antioxidant activity. Is not associated with high density lipoprotein. Prevents LDL lipid peroxidation, reverses the oxidation of mildly oxidized LDL, and inhibits the ability of MM-LDL to induce monocyte chemotaxis.


Research Area

Neuroscience

Image Data



Western blot analysis of PON2 in K562, HeLa lysates using PON2 antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using PON2 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.