

Product Name: PRPF4 Rabbit Monoclonal Antibody
Catalog #: AMRe02485



Summary

Production Name	PRPF4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	PRPF4
Alternative Names	PRP4; RP70; HPRP4; Prp4p; HPRP4P; SNRNP60
Gene ID	9128
SwissProt ID	O43172.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 58 kDa; Observed MW: 58 kDa

Background

Participates in pre-mRNA splicing. Part of the U4/U5/U6 tri-snRNP complex, one of the building blocks of the spliceosome.


Product Name: PRPF4 Rabbit Monoclonal Antibody
Catalog #: AMRe02485



Research Area

Cell Biology

Image Data



Western blot analysis of PRPF4 in K562, 3T3, HeLa lysates using PRPF4 antibody.

Note

For research use only.