

Product Name: Sall4 Rabbit Monoclonal Antibody
Catalog #: AMRe02570



Summary

Production Name	Sall4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	SALL4
Alternative Names	Zinc finger protein 797; Zinc finger protein SALL4
Gene ID	57167
SwissProt ID	Q9UJQ4.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 112 kDa; Observed MW: 112 kDa

Background

Transcription factor with a key role in the maintenance and self-renewal of embryonic and hematopoietic stem cells.


Product Name: Sall4 Rabbit Monoclonal Antibody
Catalog #: AMRe02570



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Sall4 in Daudi lysates using Sall4 antibody.

Note

For research use only.