

Summary

Production Name	SFRP4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	SFRP4
Alternative Names	FRP4; FrpHE; PYL; Sfrp4
Gene ID	6424
SwissProt ID	Q6FHJ7.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 53-55 kDa

Background

Product Name: SFRP4 Rabbit Monoclonal Antibody
Catalog #: AMRe02594




Soluble frizzled-related proteins (sFRPS) function as modulators of Wnt signaling through direct interaction with Wnts. They have a role in regulating cell growth and differentiation in specific cell types.


Research Area

Stem Cells

Image Data



Western blot analysis of SFRP4 in rat kidney, rat heart lysates using SFRP4 antibody.



Western blot analysis of SFRP4 in K562, 3T3, HeLa lysates using SFRP4 antibody.

Note

For research use only.