

**Product Name: T Box Protein 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02671**

---

## Summary

<b>Production Name</b>	T Box Protein 3 Rabbit Monoclonal Antibody
<b>Description</b>	Rabbit Monoclonal antibody
<b>Host</b>	Rabbit
<b>Application</b>	WB,IP
<b>Reactivity</b>	Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	TBX3
<b>Alternative Names</b>	T box 3; T-box protein 3; Tbx3; TBX3 ISO; UMS; XHL
<b>Gene ID</b>	6926
<b>SwissProt ID</b>	O15119.

## Application

<b>Dilution Ratio</b>	WB: 1:500-1:1000 IP: 1:20
<b>Molecular Weight</b>	Calculated MW: 79 kDa; Observed MW: 79 kDa

## Background

---

**Product Name: T Box Protein 3 Rabbit Monoclonal Antibody**  
**Catalog #: AMRe02671**

---




Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.

## Research Area

Epigenetics and Nuclear Signaling

## Image Data



Western blot analysis of Tbx3 in 3T3 lysates using T Box Protein 3 antibody.

## Note

For research use only.