

Product Name: ZNF160 Rabbit Monoclonal Antibody
Catalog #: AMRe02789



Summary

Production Name	ZNF160 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	ZNF160
Alternative Names	Zinc finger protein 160; Zinc finger protein HZF5; Zinc finger protein Kr18
Gene ID	90338
SwissProt ID	Q9HCG1.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 94 kDa; Observed MW: 74 kDa

Background

May be involved in transcriptional regulation.


Product Name: ZNF160 Rabbit Monoclonal Antibody
Catalog #: AMRe02789



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZNF160 in mouse liver lysates using ZNF160 antibody.

Note

For research use only.