

Product Name: ETS1 Rabbit Monoclonal Antibody
Catalog #: AMRe02807



Summary

Production Name	ETS1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human,Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	ETS1
Alternative Names	ETS1; EWSR2; Protein C-ets-1; p54
Gene ID	2113
SwissProt ID	P14921.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 50 kDa; Observed MW: 50 kDa

Background

Transcription factor.


Product Name: ETS1 Rabbit Monoclonal Antibody
Catalog #: AMRe02807



Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ETS1 in HeLa, A549, HL-60, U2OS lysates using ETS1 antibody.

Note

For research use only.