

Product Name: Phospho-Nucleophosmin (Thr199)
Rabbit Monoclonal Antibody
Catalog #: AMRe02863

Summary

Production Name	Phospho-Nucleophosmin (Thr199) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Rat

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	NPM1
Alternative Names	B23; NPM
Gene ID	4869
SwissProt ID	P06748.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 33 kDa; Observed MW: 38 kDa

Product Name: Phospho-Nucleophosmin (Thr199)
Rabbit Monoclonal Antibody
Catalog #: AMRe02863




Background

Involved in diverse cellular processes such as ribosome biogenesis, centrosome duplication, protein chaperoning, histone assembly, cell proliferation, and regulation of tumor suppressors p53/TP53 and ARF. Binds ribosome presumably to drive ribosome nuclear export. Associated with nucleolar ribonucleoprotein structures and bind single-stranded nucleic acids. Acts as a chaperonin for the core histones H3, H2B and H4. Stimulates APEX1 endonuclease activity on apurinic/apyrimidinic (AP) double-stranded DNA but inhibits APEX1 endonuclease activity on AP single-stranded RNA. May exert a control of APEX1 endonuclease activity within nucleoli devoted to repair AP on rDNA and the removal of oxidized rRNA molecules. In concert with BRCA2, regulates centrosome duplication. Regulates centriole duplication: phosphorylation by PLK2 is able to trigger centriole replication. Negatively regulates the activation of EIF2AK2/PKR and suppresses apoptosis through inhibition of EIF2AK2/PKR autophosphorylation. Antagonizes the inhibitory effect of ATF5 on cell proliferation and relieves ATF5-induced G2/M blockade (PubMed:22528486). In complex with MYC enhances the transcription of MYC target genes (PubMed:25956029).

Research Area


Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of Nucleophosmin (Phospho-Thr199) in C6, Ramos lysates using Phospho-Nucleophosmin (Thr199) antibody.

Product Name: Phospho-Nucleophosmin (Thr199)
Rabbit Monoclonal Antibody
Catalog #: AMRe02863



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Nucleophosmin (Phospho-Thr199) antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.