

Summary

Production Name	CAB39 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	CAB39
Alternative Names	MO25; CGI-66
Gene ID	51719
SwissProt ID	Q9Y376.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 40 kDa; Observed MW: 40 kDa

Background

Product Name: CAB39 Rabbit Monoclonal Antibody
Catalog #: AMRe02920




Component of a complex that binds and activates STK11/LKB1.


Research Area

Signal Transduction


Image Data



Western blot analysis of CAB39/MO25 in mouse kidney, rat lung lysates using CAB39 antibody.



Western blot analysis of CAB39/MO25 in Hela, A549, HL-60, U2OS, U87-MG lysates using CAB39/MO25 antibody.



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using CAB39/MO25 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: CAB39 Rabbit Monoclonal Antibody
Catalog #: AMRe02920



Note

For research use only.