

Product Name: HJURP Rabbit Monoclonal Antibody
Catalog #: AMRe03030



Summary

Production Name	HJURP Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	HJURP
Alternative Names	FAKTS; URLC9; hFLEG1
Gene ID	55355
SwissProt ID	Q8NCD3.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 84 kDa

Background

Product Name: HJURP Rabbit Monoclonal Antibody
Catalog #: AMRe03030




Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres.


Research Area

-

Image Data



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.