

Product Name: IKB epsilon Rabbit Monoclonal Antibody
Catalog #: AMRe03044



Summary

Production Name	IKB epsilon Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	NFKBIE
Alternative Names	NFKBIE; IKBE; NF-kappa-B inhibitor epsilon; NF-kappa-BIE; I-kappa-B-epsilon; IkB-E; IkB-epsilon; IkappaBepsilon
Gene ID	4794
SwissProt ID	O00221.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 53 kDa; Observed MW: 48 kDa

Background

Product Name: IKB epsilon Rabbit Monoclonal Antibody
Catalog #: AMRe03044




Inhibits NF-kappa-B by complexing with and trapping it in the cytoplasm. Inhibits DNA-binding of NF-kappa-B p50-p65 and p50-c-Rel complexes.

Research Area

Signal Transduction

Image Data



Western blot analysis of IKB epsilon in HeLa, Ramos lysates using IKB epsilon antibody.

Note

For research use only.