

Product Name: Meckelin Rabbit Monoclonal Antibody
Catalog #: AMRe03081



Summary

Production Name	Meckelin Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	TMEM67
Alternative Names	MKS3; JBTS6; NPHP11; TNEM67; MECKELIN
Gene ID	91147
SwissProt ID	Q5HYA8.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 112 kDa; Observed MW: 112 kDa

Background

Product Name: Meckelin Rabbit Monoclonal Antibody
Catalog #: AMRe03081




Required for ciliary structure and function.


Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Meckelin in HeLa, A549, HL-60, U251, U87-MG lysates using Meckelin antibody.



Immunocytochemistry analysis of Meckelin (green) in MCF-7 using Meckelin antibody, and DAPI (blue)

Note

For research use only.