

Product Name: MVK Rabbit Monoclonal antibody
Catalog #: AMRe03101



Summary

Production Name	MVK Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	MVK
Alternative Names	MVK; Mevalonate kinase; MK
Gene ID	4598
SwissProt ID	Q03426.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 42 kDa; Observed MW: 42 kDa

Background

Catalyzes the phosphorylation of mevalonate to mevalonate 5-phosphate, a key step in isoprenoid and cholesterol

Product Name: MVK Rabbit Monoclonal antibody
Catalog #: AMRe03101




biosynthesis.

Research Area

-

Image Data



Western blot analysis of MVK in HeLa, A549, HL-60, U251, U87-MG lysates using MVK antibody.

Note

For research use only.