

Product Name: NUP155 Rabbit Monoclonal Antibody
Catalog #: AMRe03123



Summary

Production Name	NUP155 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	NUP155
Alternative Names	N155; ATFB15
Gene ID	9631
SwissProt ID	O75694.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 155 kDa; Observed MW: 155 kDa

Background

Product Name: NUP155 Rabbit Monoclonal Antibody
Catalog #: AMRe03123




Essential component of nuclear pore complex.


Research Area

Tags & Cell Markers

Image Data



Immunocytochemistry analysis of NUP155 (green) in HeLa using NUP155 antibody, and DAPI (blue).



Western blot analysis of NUP155 in HeLa, A549, HL-60 lysates using NUP155 antibody.

Note

For research use only.