

Product Name: P2Y12 Rabbit Monoclonal antibody
Catalog #: AMRe03133



Summary

Production Name	P2Y12 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	P2RY12
Alternative Names	HORK3; P2Y12; ADPG-R; BDPLT8; SP1999; P2T(AC); P2Y(AC); P2Y(12)R; P2Y(ADP); P2Y(cyc)
Gene ID	64805
SwissProt ID	Q9H244.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 39 kDa; Observed MW: 52 kDa

Background

Product Name: P2Y12 Rabbit Monoclonal antibody
Catalog #: AMRe03133




Receptor for ADP and ATP coupled to G-proteins that inhibit the adenylyl cyclase second messenger system.

Research Area

Neuroscience

Image Data



Western blot analysis of P2Y12 in A549, HL-60, U251, U87-MG lysates using P2Y12 antibody.

Note

For research use only.