

Product Name: PPCEL Rabbit Monoclonal antibody
Catalog #: AMRe03158



Summary

Production Name	PPCEL Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-P,IP
Reactivity	Human,Mouse,Rat,Hamster

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	PREPL
Alternative Names	CMS22; Prolylendopeptidase-like
Gene ID	9581
SwissProt ID	Q4J6C6.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IP: 1:20
Molecular Weight	Calculated MW: 84 kDa; Observed MW: 72 kDa

Background

Product Name: PPCEL Rabbit Monoclonal antibody
Catalog #: AMRe03158




Serine peptidase whose precise substrate specificity remains unclear.


Research Area

-


Image Data



Western blot analysis of PPCEL/PREPL in mouse brain lysates using PPCEL antibody.



Western blot analysis of PPCEL/PREPL in HeLa, CHO-K1, C6 lysates using PPCEL/PREPL antibody.



Immunohistochemistry analysis of paraffin-embedded Human breast cancer using PPCEL/PREPL antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Product Name: PPCEL Rabbit Monoclonal antibody
Catalog #: AMRe03158



Note

For research use only.