

Product Name: RBM3 Rabbit Monoclonal antibody
Catalog #: AMRe03183



Summary

Production Name	RBM3 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	RBM3
Alternative Names	RNPL; IS1-RNPL
Gene ID	5935
SwissProt ID	P98179.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 17 kDa; Observed MW: 17 kDa

Background

Product Name: RBM3 Rabbit Monoclonal antibody
Catalog #: AMRe03183




Cold-inducible mRNA binding protein that enhances global protein synthesis at both physiological and mild hypothermic temperatures.


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of RBM3 in 3T3, mouse brain lysates using RBM3 antibody.



Western blot analysis of RBM3 in A549, HL-60, U251 lysates using RBM3 antibody.



Immunocytochemistry analysis of RBM3 (green) in hela using RBM3 antibody, and DAPI (blue)

Product Name: RBM3 Rabbit Monoclonal antibody
Catalog #: AMRe03183



Note

For research use only.