

Product Name: Transcription Initiation Factor IIB Rabbit Monoclonal Antibody
Catalog #: AMRe03228

Summary

Production Name	Transcription Initiation Factor IIB Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	GTF2B
Alternative Names	General transcription factor IIB; General transcription factor TFIIB; gtf2b; RNA polymerase II transcription factor IIB; S300 II; S300-II; TF IIB; TF2B; TF2B_HUMAN; TFIIB; Transcription initiation factor IIB.
Gene ID	2959
SwissProt ID	Q00403.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 35 kDa; Observed MW: 35 kDa

**Product Name: Transcription Initiation Factor IIB Rabbit
Monoclonal Antibody**
Catalog #: AMRe03228




Background

General factor that plays a major role in the activation of eukaryotic genes transcribed by RNA polymerase II.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of TFIIIB in 3T3, Raw264.7 lysates using Transcription Initiation Factor IIB antibody.



Western blot analysis of TFIIIB in HeLa, A549, HL-60, U2OS, C6 lysates using TFIIIB antibody.

Note

For research use only.