

Summary

Production Name	USP22 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purification

Immunogen

Gene Name	USP22
Alternative Names	Ubiquitin carboxyl terminal hydrolase 22; Usp22; USP3L
Gene ID	23326
SwissProt ID	Q9UPT9.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Background

Product Name: USP22 Rabbit Monoclonal Antibody
Catalog #: AMRe03253




Histone deubiquitinating component of the transcription regulatory histone acetylation (HAT) complex SAGA. Catalyzes the deubiquitination of both histones H2A and H2B, thereby acting as a coactivator.


Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of USP22 in F9 lysates using USP22 antibody.



Western blot analysis of USP22 in A549, HL-60, U2OS, C6 lysates using USP22 antibody.

Note

For research use only.