

Product Name: ZBTB7A Rabbit Monoclonal antibody
Catalog #: AMRe03260



Summary

Production Name	ZBTB7A Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ZBTB7A
Alternative Names	FBI1; LRF; Pokemon; TIP21; ZBTB7; ZBTB7A; ZNF857A
Gene ID	51341
SwissProt ID	O95365.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 61 kDa; Observed MW: 61 kDa

Background

Product Name: ZBTB7A Rabbit Monoclonal antibody
Catalog #: AMRe03260




Plays a key role in the instruction of early lymphoid progenitors to develop into B lineage by repressing T-cell instructive Notch signals (By similarity).


Research Area

Epigenetics and Nuclear Signaling


Image Data



Western blot analysis of ZBTB7A in HeLa, Ramos lysates using ZBTB7A antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using ZBTB7A antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Immunocytochemistry analysis of ZBTB7A (green) in MCF-7 using ZBTB7A antibody, and DAPI (blue)

Product Name: ZBTB7A Rabbit Monoclonal antibody
Catalog #: AMRe03260



Note

For research use only.