

Product Name: ZNF217 Rabbit Monoclonal antibody
Catalog #: AMRe03263



Summary

Production Name	ZNF217 Rabbit Monoclonal antibody
Description	Recombinant Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal Antibody
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40% Glycerol, 0.01% Sodium azide and 0.05% BSA
Purification	Affinity Purified

Immunogen

Gene Name	ZNF217
Alternative Names	ZABC1; ZN217; ZNF217
Gene ID	7764
SwissProt ID	O75362.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW: 115 kDa; Observed MW: 115 kDa

Background

Binds to the promoters of target genes and functions as repressor. Promotes cell proliferation and antagonizes cell death.

Product Name: ZNF217 Rabbit Monoclonal antibody
Catalog #: AMRe03263




Promotes phosphorylation of AKT1 at 'Ser-473'.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of ZNF217 in A549, HL-60, U251 lysates using ZNF217 antibody.

Note

For research use only.