

Product Name: PGC1 beta Rabbit Monoclonal Antibody
Catalog #: AMRe03945



Summary

Production Name	PGC1 beta Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IP
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	PPARGC1B
Alternative Names	PERC; PGC1; PPARGC1; Ppargc1b
Gene ID	133522
SwissProt ID	Q86YN6.

Application

Dilution Ratio	WB: 1:500-1:1000 IP: 1:20
Molecular Weight	Calculated MW:113 kDa;Observed MW: 113 kDa

Background

Plays a role of stimulator of transcription factors and nuclear receptors activities. Activates transcriptional activity of estrogen

Product Name: PGC1 beta Rabbit Monoclonal Antibody
Catalog #: AMRe03945




receptor alpha, nuclear respiratory factor 1 (NRF1) and glucocorticoid receptor in the presence of glucocorticoids.

Research Area

Signal Transduction

Image Data



Western blot analysis of PGC1 beta in HeLa, MCF-7 lysates using PGC1 beta antibody.

Note

For research use only.