

Product Name: PPIC Rabbit Monoclonal Antibody
Catalog #: AMRe03949



Summary

Production Name	PPIC Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	PPIC
Alternative Names	PPIC,CYPC,Rotamase C,Cyclophilin C,Peptidyl-prolyl cis-trans isomerase C
Gene ID	5480.0
SwissProt ID	P45877.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW:23 kDa;Observed MW: 23 kDa

Background

PPIase that catalyzes the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and may therefore assist

Product Name: PPIC Rabbit Monoclonal Antibody
Catalog #: AMRe03949




protein folding.

Research Area

Immunology

Image Data



Western blot analysis of PPIC in HCT116 lysates using PPIC antibody.

Note

For research use only.