

Product Name: CAD Rabbit Monoclonal Antibody
Catalog #: AMRe03951



Summary

Production Name	CAD Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	CAD
Alternative Names	CAD; CAD protein
Gene ID	790
SwissProt ID	P27708.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW:243 kDa;Observed MW: 243 kDa

Background

Carbamoyl phosphate synthetase-aspartate carbamoyltransferase-dihydroorotase (CAD) is a multifunctional protein that

Product Name: CAD Rabbit Monoclonal Antibody
Catalog #: AMRe03951




initiates and regulates mammalian de novo pyrimidine biosynthesis.

Research Area

Cell Biology

Image Data



Western blot analysis of CAD in 293 lysates using CAD antibody.

Note

For research use only.