

Product Name: Yes1 Rabbit Monoclonal Antibody
Catalog #: AMRe03958



Summary

Production Name	Yes1 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	YES1
Alternative Names	YES1; YES; Tyrosine-protein kinase Yes; Proto-oncogene c-Yes; p61-Yes
Gene ID	7525
SwissProt ID	P07947.

Application

Dilution Ratio	WB: 1:500-1:1000
Molecular Weight	Calculated MW:61 kDa;Observed MW: 61 kDa

Background

Promotes infectivity of Neisseria gonorrhoeae in epithelial cells by phosphorylating MCP/CD46.


Product Name: Yes1 Rabbit Monoclonal Antibody
Catalog #: AMRe03958



Research Area

Signal Transduction

Image Data



Western blot analysis of Yes1 in SW480, Jurkat lysates using Yes1 antibody.

Note

For research use only.