

Product Name: MRPS35 Rabbit Monoclonal Antibody
Catalog #: AMRe03968



Summary

Production Name	MRPS35 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	MRPS35
Alternative Names	mitochondrial ribosomal protein S35; MDS023; MRPS28; MRP-S28; HDCMD11P
Gene ID	60488
SwissProt ID	P82673.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100
Molecular Weight	Calculated MW:37 kDa;Observed MW: 37 kDa

Background

Mammalian mitochondrial ribosomal proteins are encoded by nuclear genes and help in protein synthesis within the

Product Name: MRPS35 Rabbit Monoclonal Antibody
Catalog #: AMRe03968




mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit.

Research Area

-

Image Data



Western blot analysis of MRPS35 in HeLa ,Jurkat lysates using MRPS35 antibody.

Note

For research use only.