

Product Name: Pirh2 Rabbit Monoclonal Antibody
Catalog #: AMRe03971



Summary

Production Name	Pirh2 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	RCHY1
Alternative Names	RCHY1; ARNIP; CHIMP; PIRH2; RNF199; ZNF363; RING finger and CHY zinc finger domain-containing protein 1; Androgen receptor N-terminal-interacting protein; CH-rich-interacting match with PLAG1; E3 ubiquitin-protein ligase Pirh2; RING finger
Gene ID	25898
SwissProt ID	Q96PM5.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW:30 kDa;Observed MW: 30 kDa

Product Name: Pirh2 Rabbit Monoclonal Antibody
Catalog #: AMRe03971




Background

Mediates E3-dependent ubiquitination and proteasomal degradation of target proteins, including p53/TP53, HDAC1 and CDKN1B. Preferentially acts on tetrameric p53/TP53.

Research Area

Cell Biology

Image Data



Western blot analysis of Pirh2 in HeLa lysates using Pirh2 antibody.

Note

For research use only.