

Product Name: RPLP0 Rabbit Monoclonal Antibody
Catalog #: AMRe03976



Summary

Production Name	RPLP0 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,IHC-F,IHC-P,ICC/IF
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	RPLP0
Alternative Names	Ribosomal protein; large; P0; P0; LP0; L10E; RPP0; PRLP0
Gene ID	6175
SwissProt ID	P05388.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200
Molecular Weight	Calculated MW:34 kDa;Observed MW: 34 kDa

Background

Ribosomal protein P0 is the functional equivalent of E.


Product Name: RPLP0 Rabbit Monoclonal Antibody
Catalog #: AMRe03976



Research Area

-

Image Data



Western blot analysis of RPLP0 in 293, HeLa lysates using RPLP0 antibody.

Note

For research use only.