

Product Name: Ferritin Light Chain Rabbit Monoclonal Antibody
Catalog #: AMRe04005



Summary

Production Name	Ferritin Light Chain Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB,FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	FTL
Alternative Names	Ferritin L chain; Ferritin L subunit; Ferritin light chain; Ferritin light polypeptide; ferritin light polypeptide like 3; FTL; LFTD; NBIA 3;
Gene ID	2512
SwissProt ID	P02792.

Application

Dilution Ratio	WB: 1:500-1:1000 FC: 1:50-1:100
Molecular Weight	Calculated MW:20 kDa;Observed MW: 20 kDa

Product Name: Ferritin Light Chain Rabbit Monoclonal Antibody
Catalog #: AMRe04005




Background

Stores iron in a soluble, non-toxic, readily available form. Important for iron homeostasis. Iron is taken up in the ferrous form and deposited as ferric hydroxides after oxidation. Also plays a role in delivery of iron to cells. Mediates iron uptake in capsule cells of the developing kidney.

Research Area

Neuroscience

Image Data



Western blot analysis of Ferritin Light Chain in HepG2, HeLa lysates using Ferritin Light Chain antibody.

Note

For research use only.