

**Product Name: MonoMethyl-Histone H3 (Arg2) Rabbit
Monoclonal Antibody
Catalog #: AMRe04012**



Summary

Production Name	MonoMethyl-Histone H3 (Arg2) Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Monomethylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	H3C1
Alternative Names	H3R2me; H3/j; H3C1; H3C2; H3C3; H3C4; H3C6; H3C7; H3C8; H3FJ; H3C10; H3C11; HIST1H3J
Gene ID	8350
SwissProt ID	P68431.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW:15 kDa;Observed MW: 17 kDa

**Product Name: MonoMethyl-Histone H3 (Arg2) Rabbit
Monoclonal Antibody
Catalog #: AMRe04012**




Background

H3 Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.

Research Area

Epigenetics and Nuclear Signaling

Image Data



Western blot analysis of MonoMethyl-Histone H3 (Arg2) in HeLa, 3T3 lysates using MonoMethyl-Histone H3 (Arg2) antibody.

Note

For research use only.