

Product Name: Annexin 4/ANXA4 Rabbit Monoclonal Antibody
Catalog #: AMRe04025



Summary

Production Name	Annexin 4/ANXA4 Rabbit Monoclonal Antibody
Description	Rabbit Monoclonal antibody
Host	Rabbit
Application	WB, ICC/IF, FC
Reactivity	Human, Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Liquid in 50mM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.
Purification	Affinity Purification

Immunogen

Gene Name	ANXA4
Alternative Names	AlV; AnnexinA4; AnnexinIV; ANX 4; ANX A4; ANX4; ANXA4; Chromobindin4; EndonexinI; LipocortinIV; P32.5; P33/41; PAPII; PIG28; PP4X; Protein II; Xanx-4; ZAP36
Gene ID	307
SwissProt ID	P09525.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200 FC: 1:50-1:100
Molecular Weight	Calculated MW:36 kDa;Observed MW: 36 kDa

Product Name: Annexin 4/ANXA4 Rabbit Monoclonal Antibody
Catalog #: AMRe04025




Background

Calcium/phospholipid-binding protein which promotes membrane fusion and is involved in exocytosis.

Research Area

Cardiovascular

Image Data



Western blot analysis of Annexin 4/ANXA4 in Human placenta, HepG2 lysates using Annexin 4/ANXA4 antibody.

Note

For research use only.